

CE  
2460



ENHANCED PROTECTION WITH GREATER VISIBILITY

# BIO-HYALUR<sup>TM</sup>CS

Sodium Hyaluronate 2% &  
Sodium Chondroitin Sulfate 2% Ophthalmic Solution

biotech  
VISION CARE

# BIO-HYALUR<sup>TM</sup>CS

Sodium Hyaluronate 2% &  
Sodium Chondroitin Sulfate 2% Ophthalmic Solution

## Features

- Combination of Sodium Hyaluronate + Sodium Chondroitin Sulfate
- Excellent protection to Corneal Endothelium
- Excellent Chamber maintenance & stability
- Very quiet eyes post-operatively

## Specifications

Raw Material	Sodium Hyaluronate content (2.0% W/V) & Sodium Chondroitin Sulfate content (2.0% W/V)
Viscosity	@ 1 Shear Rate - 40,000 ± 15,000 mPas
pH	6.8 to 7.6
Osmolarity	340 - 400 mOsm/kg
Shelf life	2 years from manufacture
Storage	Store between 2°C to 25°C. Do not Freeze. Protect from light
Application / Intended Use	Anterior segment procedures



biotech  
VISION CARE