

CE
2460

biotech

VISION CARE

CATARACT SURGERY SOLUTIONS

www.biotechhealthcare.com

Eye Visc

Hydroxy Propyl Methyl Cellulose (HPMC) Ophthalmic Solution USP 2.0% W/V


- Highly retentive, Dispersive viscoelastic
- Unique ability to offer good space maintenance and excellent tissue protection throughout the procedure
- Possesses a high dynamic viscoelasticity and coating ability
- Improves the depth of anterior chamber allowing easy and safe implantation of IOLs
- Excellent ease of use
- Viscosity 4500-9500 mPas

Indications:

- Eye Visc is indicated when ever protection & lubrication on the delicate cells & tissues are needed especially in ophthalmic procedure including Anterior Segment surgery, Cataract surgery, IOL implantation, Corneal transplantation & Glaucoma surgery.

Packaging:

- Eyevisc is supplied 2.0 ml in sterile PFS





biotech

VISION CARE

BIOTECH VISION CARE PVT LTD | Block 1, Abhishree Corporate Park, Opp. Swagat Bunglow BRTS Stop
Bopal - Ambli Road, Ahmedabad - 380 058 Gujarat India. | Tel: +91 79 66823000

Email: indiasales@biotechhealthcare.com | intlsales@biotechhealthcare.com | www.biotechhealthcare.com