


BIO SPEARS

Sterile Surgical Spears



TOOLS FOR FLUID CONTROL

BIO SPEARS

Sterile Surgical Spears



BIO SPEARS (Sterile PVA Eye Spears)

- Designed to manipulate tissue during LASIK and refractive procedures
- Offers doctors with superior fibre free products for unparalleled fluid management and control
- Lint free, high density PVA formulation and virtually particulate free
- Highly absorbent, single use, fast wicking and packaged sterile
- Maintains rigidity during wicking process
- Problems relating to cotton fibres - left - overs are completely avoided

